

SILICONE TRACHEOSTOMY TUBE

Silicone for Enhanced Patient Comfort and Great
Biocompatibility



**ANESTHETIC
SERIES**

SILICONE TRACHEOSTOMY TUBE

Fortune® Silicone Tracheostomy Tube Features:

- Medical grade silicone for superior
- biocompatibility.
- Curved smooth tube and mounted Obturator for easy insertion.
- Standard 15 mm connector.
- X-ray opaque line.
- All cuff models comes with a sensitive pilot balloon allows confirmation of the cuff status.
- Low pressure cuff model provides minimum pressure against trachea wall.
- Neck strap for fixation.
- 15 mm Disconnect Wedge.
- Sterilized double packaging.

- A** Tracheostomy Tube, Cuffless
- B** Tracheostomy Tube, OC Type
- C** Tracheostomy Tube, LC Type
- D** Disconnect Wedge
- E** Neck Strap




SILICONE TRACHEOSTOMY TUBE

SPECIFICATIONS



Disconnect Wedge:

REF. No.	Size	Description
1500-0000	15 mm	For disconnect 15 mm connector




Neck Strap:

REF. No.	Size	Description
1550-0001	Universal	Soft fabric Velcro



Silicone Tracheostomy Tube, Cuffless:

Product	REF. No.	Size / I.D.	O.D.	Length	Description
Tracheostomy Tube, Cuffless	1561-0050	5.0 mm	7.3 mm	60.0 mm	
	1561-0055	5.5 mm	8.0 mm	62.5 mm	
	1561-0060	6.0 mm	8.7 mm	65.0 mm	
	1561-0065	6.5 mm	9.3 mm	67.5 mm	- Cuffless
	1561-0070	7.0 mm	10.0 mm	70.0 mm	- X-ray opaque line
	1561-0075	7.5 mm	10.7 mm	72.5 mm	- Obturator
	1561-0080	8.0 mm	11.3 mm	75.0 mm	- Neck strap
	1561-0085	8.5 mm	12.0 mm	77.5 mm	- Disconnect wedge
	1561-0090	9.0 mm	12.7 mm	80.0 mm	
	1561-0095	9.5 mm	13.3 mm	82.5 mm	
	1561-0100	10.0 mm	14.0 mm	85.0 mm	




SILICONE TRACHEOSTOMY TUBE

Silicone Tracheostomy, OC Type:


Product	REF. No.	Size / I.D.	O.D.	Length	Description
Tracheostomy Tube, OC Type	1560-0050	5.0 mm	7.3 mm	60.0 mm	
	1560-0055	5.5 mm	8.0 mm	62.5 mm	
	1560-0060	6.0 mm	8.7 mm	65.0 mm	
	1560-0065	6.5 mm	9.3 mm	67.5 mm	- Ordinary Cuff
	1560-0070	7.0 mm	10.0 mm	70.0 mm	- X-ray opaque line
	1560-0075	7.5 mm	10.7 mm	72.5 mm	- Obturator
	1560-0080	8.0 mm	11.3 mm	75.0 mm	- Neck strap
	1560-0085	8.5 mm	12.0 mm	77.5 mm	- Disconnect wedge
	1560-0090	9.0 mm	12.7 mm	80.0 mm	
	1560-0095	9.5 mm	13.3 mm	82.5 mm	
	1560-0100	10.0 mm	14.0 mm	85.0 mm	

Silicone Tracheostomy, LC Type:

Product	REF. No.	Size / I.D.	O.D.	Length	Description
Tracheostomy Tube, LC Type	1565-0050	5.0 mm	7.3 mm	60.0 mm	
	1565-0055	5.5 mm	8.0 mm	62.5 mm	
	1565-0060	6.0 mm	8.7 mm	65.0 mm	
	1565-0065	6.5 mm	9.3 mm	67.5 mm	- Low Pressure Cuff
	1565-0070	7.0 mm	10.0 mm	70.0 mm	- X-ray opaque line
	1565-0075	7.5 mm	10.7 mm	72.5 mm	- Obturator
	1565-0080	8.0 mm	11.3 mm	75.0 mm	- Neck strap
	1565-0085	8.5 mm	12.0 mm	77.5 mm	- Disconnect wedge
	1565-0090	9.0 mm	12.7 mm	80.0 mm	
	1565-0095	9.5 mm	13.3 mm	82.5 mm	
	1565-0100	10.0 mm	14.0 mm	85.0 mm	



1560-OXXX



1565-OXXX



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com