







# TRACHEAL TUBE FIXATION DEVICE

## **INTENDED USE:**

The Fortune® tracheal tube fixation device (hereinafter referred to as the fixer) is intended to be used to secure a tracheal tube in an appropriate position. It is easy to operate and enhances patient comfort.

#### **FEATURES:**

- 1. Easy to operate: quick to mount; efficiently prevents biting of the tube.
- 2. The anti-ulceration pad reduces the risk of ulcerated lips.
- 3. Smaller and shorter than artificial airways, the fixer's small size reduces the risk of damage to the patient's oral cavity.
- 4. The fixer's adjustable strap reduces the risk of skin damage. The strap is designed with patient comfort in mind.




Patented design wins 2013 Taipei Int'l Invention
Show & Technomart Gold Medal Award

# **SPECIFICATIONS**

Designed to match tracheal tubes with outer diameter (O.D.) from 8.0 mm to 12.0 mm.

Fortune Tracheal tube		size	O.D. (mm)	I.D. (mm)
	0060	26FR	8.7	6.0
1520 series, 1525 series, 1530 series, 1535 series, 1540 series, 1545 series, 1550 series, 1555 series.	0065	28FR	9.3	6.5
	0070	30FR	10.0	7.0
	0075	32FR	10.7	7.5
	0080	34FR	11.3	8.0
	0085	36FR	12.0	8.5
	0090	38FR	12.7	9.0



1500-0101

## **Head Strap:**

REF. No.	Size	Description
1500-0111	Universal	Soft fabric Velcro

