

SILICONE WOUND DRAINAGE SYSTEM

Silicone for great biocompatibility
and fast recovery

- Silicone C.W.V. Reservoir...5-2
- Collection Bag...5-3
- Wound Drainage Trocar...5-4
- Silicone Flat Drain...5-5
- Silicone Round Drain...5-8
- Silicone Vac Drain...5-11
- Silicone Endo Drain...5-15
- Silicone Hubless Endo Drain...5-19
- Endo Connector...5-23

PRODUCT SET

- Wound Drainage Tube with FixPad®...5-24

SILICONE WOUND DRAINAGE SYSTEM



Post-operative wound drainage can be safe and easy.

Made of medical grade silicone, FortuneMed® provides a complete line of C.W.V. (closed wound vacuum) Reservoirs, Drains, Trocars and Collection Bags for different surgical applications.

Why Silicone?

The advantages of silicone provide superior biocompatibility, flexibility and great comfort to patient. For this reason, Fortune® silicone wound drain systems minimize tissue adhesion and trauma during both healing and removal, resulting in extremely lower tissue response compared to drains made of other elastomers, such as PVC and latex.

- A Silicone C.W.V. Reservoir**
- B Collection Bag**
- C Trocar**
- D Silicone Flat Drain**
- E Silicone Round Drain**
- F Silicone Vac Drain**
- G Silicone Endo Drain**
- H Silicone Hubless Endo Drain**
- I Endo Connector**




*Data based on 102 female patients operated, 52 using PVC wound drainage system and 50 using Fortune® Silicone Wound Drain System. Drain pull-out time for measuring fluid drawn volume is 6.18 days for PVC system and 5.18 days for silicone system.

SILICONE C.W.V. RESERVOIR

Fortune® Silicone C.W.V. Reservoir Features:

- Medical grade silicone for superior biocompatibility.
- Low pressure design provides continuous gentle suction and minimizes tissue trauma.
- Volume marking provides easy visual observation.
- Flat bottom design for easy positioning.
- Non-reflex valve.
- Hanger and bedside clip for secured attachment.
- A great variety of sizes and types are available.
- Optional extra inlet is available on 200 & 400 cc models to connect second tube.
- Bottom outlet models are capable to connect with Collection Bag.
- Sterilization by E.O. gas.


A Silicone C.W.V. Reservoir



SPECIFICATIONS

C.W.V. Reservoir:

REF. No.	Size	Inlet	Description
2013-0610	100 cc	1	
2013-0615	150 cc		Top outlet, Bedside clip
2013-0620	200 cc	1+1	
2013-0640	400 cc		
2013-1640	400 cc	1+1	Top outlet, female luer w/cap Bedside clip
2013-0910	100 cc	1	
2013-0915	150 cc		Bottom outlet w/ large bore male connector. Pinch clamp & cap. Bedside clip
2013-0920	200 cc	1+1	
2013-0940	400 cc		



COLLECTION BAG

Fortune® Collection Bag Features:

- Large bore female connector for connecting with bottom outlet Reservoir.
- Volume marking provides easy visual observation.
- The bag ties provide a secure attachment to the bedside.
- Non-reflex valve design prevents bio-hazard from medical personnel and avoids environmental pollution.


A Collection Bag




SPECIFICATIONS

Collection Bag:


REF. No.	Size	Description
2013-9920	200 cc	- Large bore female connector
2013-9940	400 cc	- Pinch clamp & cap
2013-9980	800 cc	- Bag ties



2013-9920



2013-9940




2013-9980

Bag System: Silicone C.W.V. Reservoir (Bottom outlet) + Collection Bag

REF. No.	Description
0910-9920	100 cc Reservoir, 200 cc Bag
0910-9940	100 cc Reservoir, 400 cc Bag
0915-9920	150 cc Reservoir, 200 cc Bag
0915-9940	150 cc Reservoir, 400 cc Bag
0920-9940	200 cc Reservoir, 400 cc Bag
0920-9980	200 cc Reservoir, 800 cc Bag
0940-9940	400 cc Reservoir, 400 cc Bag
0940-9980	400 cc Reservoir, 800 cc Bag

*The **Bag System** can also be supplemented with the required drain tube materials to form a **Drain Bag Set**.



0910-9920


WOUND DRAINAGE TROCAR

Trocar Features:

- Made of stainless steel.
- Sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

A Trocar

B Bendable Trocar



Bendable Trocar Features:

- Made of stainless steel.
- Bendable, sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

SPECIFICATIONS

Trocar:


REF. No.	Size	Description
2015-9008	8 FR	
2015-9010	10 FR	
2015-9016	16 FR	- Stainless steel - Protection cover
2015-9019	19 FR	
2015-9020	20 FR	
2015-9021	21 FR	

Bendable Trocar:

2015-9110	10FR	
2015-9116	16FR	- Bendable - Stainless steel
2015-9119	19FR	- Protection cover
2015-9121	21FR	



2015-9021



2015-91xx

SILICONE FLAT DRAIN

Fortune® Silicone Flat Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque effectiveness.
- Soft yet kink-resistant structure to minimize clog or obstruct.
- Clear connecting tube provides inspection and drainage visibility.
- A great variety of sizes and types are available.
- Sterilization by E.O. gas.
- Various combination is available.

- A Silicone Flat Drain**
- B Silicone Flat Drain Kit**
- C Silicone Flat Drain Set***
- D Silicone Flat Drain System***
- E Silicone Flat Drain Bag System***



SPECIFICATIONS

Flat Drain:

REF. No.	Size	Description
2014-0051	2 x 4.5 mm	
2014-0071	3 x 7 mm	3/4 perforation on 200 mm drain X-ray opaque effective
2014-0101	4 x 10 mm	
2014-0131	5 x 13 mm	
2014-0052	2 x 4.5 mm	
2014-0072	3 x 7 mm	Full perforation on 200 mm drain X-ray opaque effective
2014-0102	4 x 10 mm	
2014-0132	5 x 13 mm	

- Models with depth markings options are available.



2014-OXX1 3/4 Perforation




2014-OXX2 Full Perforation


Flat Drain Kit:

REF. No.	Size	Description
2014-9051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, attached Trocar
2014-9071	3 x 7 mm	7 mm 3/4 Flat Drain, attached Trocar
2014-9101	4 x 10 mm	10 mm 3/4 Flat Drain, attached Trocar
2014-9131	5 x 13 mm	13 mm 3/4 Flat Drain, attached Trocar
2014-9052	2 x 4.5 mm	4.5 mm full Flat Drain, attached Trocar
2014-9072	3 x 7 mm	7 mm full Flat Drain, attached Trocar
2014-9102	4 x 10 mm	10 mm full Flat Drain, attached Trocar
2014-9132	5 x 13 mm	13 mm full Flat Drain, attached Trocar

- Models with depth markings options are available.



2014-9XX1 3/4 Flat Drain with Trocar




2014-9XX2 Full Flat Drain with Trocar

SILICONE FLAT DRAIN


Flat Drain Set*:

REF. No.	Size	Description
0610-0051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0071	3 x 7 mm	7 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0101	4 x 10 mm	10 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0131	5 x 13 mm	13 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0052	2 x 4.5 mm	4.5 mm full Flat Drain, 100 cc Reservoir
● 0610-0072	3 x 7 mm	7 mm full Flat Drain, 100 cc Reservoir
● 0610-0102	4 x 10 mm	10 mm full Flat Drain, 100 cc Reservoir
● 0610-0132	5 x 13 mm	13 mm full Flat Drain, 100 cc Reservoir

● Models with depth markings options are available.




0610-OXX1 3/4 Perforation Flat Drain Set




0610-OXX2 Full Perforation Flat Drain Set

Flat Drain System*:

REF. No.	Size	Description
0610-9051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9071	3 x 7 mm	7 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9101	4 x 10 mm	10 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9131	5 x 13 mm	13 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9052	2 x 4.5 mm	4.5 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9072	3 x 7 mm	7 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9102	4 x 10 mm	10 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9132	5 x 13 mm	13 mm full Flat Drain, Trocar, 100 cc Reservoir



0610-9XX1 3/4 Flat Drain System



0610-9XX2 Full Flat Drain System


SILICONE FLAT DRAIN

Flat Drain Bag System:


REF. No.	Size	Description
0910-905120	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-907120	3 x 7 mm	7 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-910120	4 x 10 mm	10 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-913120	5 x 13 mm	13 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-905220	2 x 4.5 mm	4.5 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-907220	3 x 7 mm	7 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-910220	4 x 10 mm	10 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-913220	5 x 13 mm	13 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag

* Flat Drain Set, Flat Drain System and Flat Drain Bag System are available with 150 cc, 200 cc or 400 cc Reservoir. For detailed model numbers, please see the attached appendix.

* Flat Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-9X0120 3/4 Flat Drain Bag System




0910-9X0220 Full Flat Drain Bag System

SILICONE ROUND DRAIN

Fortune® Silicone Round Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque line.
- Rigid and flexible structure for kinking resistance.
- Clear tubing provides inspection and drainage visibility.
- A great variety of sizes and types are available.
- Sterilization by E.O. gas.
- Various combination is available.

- A Silicone Round Drain**
- B Silicone Round Drain Kit**
- C Silicone Round Drain Set***
- D Silicone Round Drain System***
- E Silicone Round Drain Bag System***



SPECIFICATIONS

Round Drain:

REF. No.	Size	Description
2015-0208	8 FR	
2015-0210	10 FR	
2015-0212	12 FR	- End perforation
2015-0214	14 FR	- 900 mm,
2015-0216	16 FR	- X-ray opaque line
2015-0218	18 FR	
2015-0220	20 FR	
2015-0308	8 FR	
2015-0310	10 FR	
2015-0312	12 FR	- Middle perforation
2015-0314	14 FR	- 1100 mm
2015-0316	16 FR	- X-ray opaque line
2015-0318	18 FR	
2015-0320	20 FR	




2015-02XX Silicone Round Drain




2015-03XX Silicone Round Drain

Round Drain Kit:

REF. No.	Size	Description
2015-9208	8 FR	8 FR end perf. Round, attached Trocar
2015-9210	10 FR	10 FR end perf. Round, attached Trocar
2015-9214	14 FR	14 FR end perf. Round, attached Trocar
2015-9216	16 FR	16 FR end perf. Round, attached Trocar
2015-9218	18 FR	18 FR end perf. Round, attached Trocar
2015-9220	20 FR	20 FR end perf. Round, attached Trocar
2015-9308	8 FR	8 FR mid. perf. Round, attached Trocar
2015-9310	10 FR	10 FR mid. perf. Round, attached Trocar
2015-9314	14 FR	14 FR mid. perf. Round, attached Trocar
2015-9316	16 FR	16 FR mid. perf. Round, attached Trocar
2015-9318	18 FR	18 FR mid. perf. Round, attached Trocar
2015-9320	20 FR	20 FR mid. perf. Round, attached Trocar



2015-92XX
Silicone Round Drain Kit




2015-93XX
Silicone Round Drain Kit


SILICONE ROUND DRAIN

Round Drain Set*:

REF. No.	Size	Description
0610-0208	8 FR	8 FR end perf. Round, 100 cc Reservoir
0610-0210	10 FR	10 FR end perf. Round, 100 cc Reservoir
0610-0212	12 FR	12 FR end perf. Round, 100 cc Reservoir
0610-0214	14 FR	14 FR end perf. Round, 100 cc Reservoir
0610-0216	16 FR	16 FR end perf. Round, 100 cc Reservoir
0610-0218	18 FR	18 FR end perf. Round, 100 cc Reservoir
0610-0220	20 FR	20 FR end perf. Round, 100 cc Reservoir
0610-0308	8 FR	8 FR mid. perf. Round, 100 cc Reservoir
0610-0310	10 FR	10 FR mid. perf. Round, 100 cc Reservoir
0610-0312	12 FR	12 FR mid. perf. Round, 100 cc Reservoir
0610-0314	14 FR	14 FR mid. perf. Round, 100 cc Reservoir
0610-0316	16 FR	16 FR mid. perf. Round, 100 cc Reservoir
0610-0318	18 FR	18 FR mid. perf. Round, 100 cc Reservoir
0610-0320	20 FR	20 FR mid. perf. Round, 100 cc Reservoir




0610-02XX Silicone Round Drain Set




0610-03XX Silicone Round Drain Set

Round Drain System*:

REF. No.	Size	Description
0610-9208	8 FR	8 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9210	10 FR	10 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9214	14 FR	14 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9216	16 FR	16 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9218	18 FR	18 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9220	20 FR	20 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9308	8 FR	8 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9310	10 FR	10 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9314	14 FR	14 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9316	16 FR	16 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9318	18 FR	18 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9320	20 FR	20 FR mid. perf. Round, Trocar, 100 cc Reservoir



0610-92XX Silicone Round Drain System



0610-93XX Silicone Round Drain System


SILICONE ROUND DRAIN

Round Drain Bag System*:


REF. No.	Size	Description
0910-920820	8 FR	8 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921020	10 FR	10 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921420	14 FR	14 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921620	16 FR	16 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921820	18 FR	18 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-922020	20 FR	20 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-930820	8 FR	8 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931020	10 FR	10 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931420	14 FR	14 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931620	16 FR	16 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931820	18 FR	18 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-932020	20 FR	20 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag

* Round Drain Set, Round Drain System and Round Drain Bag System are available with 150 cc, 200 cc or 400 cc Reservoir. For detailed model numbers, please see the attached appendix.

* Round Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-92XX20 Silicone Round Drain Bag System





0910-93XX20 Silicone Round Drain Bag System

SILICONE VAC DRAIN

Fortune® Silicone Vac Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque effectiveness.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- A great variety of sizes and types are available.
- Sterilization by E.O. gas.
- Various combination is available.

- A** Vac Drain
- B** Vac Drain Kit
- C** Vac Drain Set*
- D** Vac Drain System*
- E** Vac Drain Bag System*




SILICONE VAC DRAIN


SPECIFICATIONS

Vac Drain:

REF. No.	Size	Description
2014-0410	10 FR	- Full round Vac Drain
2014-0415	15 FR	- 200 mm round duct drain, 600 mm round tube
2014-0419	19 FR	- X-ray opaque effective
2014-0421	21 FR	- Full flat Vac Drain
2014-0408	3.5 x 8 mm	- 200 mm flat duct drain, 600 mm round tube
2014-0411	4 x 11 mm	- X-ray opaque effective




2014-0410/15/19/21




2014-0408/11

Vac Drain Kit:

REF. No.	Size	Description
2014-9410	10 FR	10 FR full Vac, attached Trocar
2014-9415	15 FR	15 FR full Vac, attached Trocar
2014-9419	19 FR	19 FR full Vac, attached Trocar
2014-9421	21 FR	21 FR full Vac, attached Trocar
2014-9408	3.5 x 8 mm	8 mm full Vac, attached Trocar
2014-9411	4 x 11 mm	11 mm full Vac, attached Trocar



2014-9410/15/19/21




2014-9408/11


SILICONE VAC DRAIN

Vac Drain Set*:

REF. No.	Size	Description
0610-0410	10 FR	10 FR full Vac, 100 cc Reservoir
0610-0415	15 FR	15 FR full Vac, 100 cc Reservoir
0610-0419	19 FR	19 FR full Vac, 100 cc Reservoir
0610-0421	21 FR	21 FR full Vac, 100 cc Reservoir
0610-0408	3.5 x 8 mm	8 mm full Vac, 100 cc Reservoir
0610-0411	4 x 11 mm	11 mm full Vac, 100 cc Reservoir




0610-0410/15/19/21




0610-0408/11

Vac Drain System*:

REF. No.	Size	Description
0610-9410	10 FR	10 FR full Vac, Trocar, 100 cc Reservoir
0610-9415	15 FR	15 FR full Vac, Trocar, 100 cc Reservoir
0610-9419	19 FR	19 FR full Vac, Trocar, 100 cc Reservoir
0610-9421	21 FR	21 FR full Vac, Trocar, 100 cc Reservoir
0610-9408	3.5 x 8 mm	8 mm full Vac, Trocar, 100 cc Reservoir
0610-9411	4 x 11 mm	11 mm full Vac, Trocar, 100 cc Reservoir



0610-9410/15/19/21



0610-9408/11

GENERAL SURGICAL SERIES


SILICONE VAC DRAIN

Vac Drain Bag System*:

REF. No.	Size	Description
0910-941020	10 FR	10 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941520	15 FR	15 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941920	19 FR	19 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-942120	21 FR	21 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-940820	3.5 x 8 mm	8 mm full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941120	4 x 11 mm	11 mm full Vac, Trocar, 100 cc Reservoir, 200 cc Bag

* Vac Drain Set, Vac Drain System and Vac Drain Bag System are available with 150 cc, 200 cc or 400 cc Reservoir. For detailed model numbers, please see the attached appendix.

* Vac Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



SILICONE ENDO DRAIN

Fortune® Silicone Endo Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque line.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- A great variety of sizes is available.
- Sterilization by E.O. gas.
- Various combination is available.

A Endo Drain

B Endo Drain Kit

C Endo Drain Set*

D Endo Drain System*


E Endo Drain Bag System*



SPECIFICATIONS

Endo Drain:


REF. No.	Size	Description
2015-0610	10 FR	
2015-0615	15 FR	
2015-0619	19 FR	- Transparent - 1200 mm
2015-0621	21 FR	- X-ray opaque line
2015-0624	24 FR	
2015-0629	29 FR	



2015-06XX Silicone Endo Drain

Endo Drain Kit:

REF. No.	Size	Description
2015-9610	10 FR	10 FR Endo, attached Trocar
2015-9615	15 FR	15 FR Endo, attached Trocar
2015-9619	19 FR	19 FR Endo, attached Trocar
2015-9621	21 FR	21 FR Endo, attached Trocar




2015-96XX Silicone Endo Drain Kit

GENERAL SURGICAL SERIES

SILICONE ENDO DRAIN


Endo Drain Set*

REF. No.	Size	Description
0610-0610	10 FR	10 FR Endo, 100 cc Reservoir
0610-0615	15 FR	15 FR Endo, 100 cc Reservoir
0610-0619	19 FR	19 FR Endo, 100 cc Reservoir
0610-0621	21 FR	21 FR Endo, 100 cc Reservoir



Endo Drain System*

REF. No.	Size	Description
0610-9610	10 FR	10 FR Endo, Trocar, 100 cc Reservoir
0610-9615	15 FR	15 FR Endo, Trocar, 100 cc Reservoir
0610-9619	19 FR	19 FR Endo, Trocar, 100 cc Reservoir
0610-9621	21 FR	21 FR Endo, Trocar, 100 cc Reservoir




Drain Bag System*

REF. No.	Size	Description
0910-961020	10 FR	10 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-961520	15 FR	15 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-961920	19 FR	19 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-962120	21 FR	21 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag

* Endo Drain Set, Endo Drain System and Endo Drain Bag System are available with 150 cc, 200 cc or 400 cc Reservoir. For detailed model numbers, please see the attached appendix.

* Endo Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



SILICONE ENDO DRAIN (TROCAR/BENDABLE TROCAR + J CONNECTOR)

Trocar Features: Please refer to the introduction on P5-4.

J Connector Features:


J Connector is designed to connect Endo Drain with the size of 10FR/15FR/19FR/21FR and Innerspring Reservoir. It is made of PVC which is a Non-Toxic material, with a little bendable character to reduce the risk of fracture. The side of J Connector designed to connect with Endo Drain is tapering, making assembly and disassembly more easily.



SPECIFICATIONS

Endo Drain Kit (Trocar + J Connector):

2015-9610J	10FR	10 FR Endo, attached Trocar, J Connector
2015-9615J	15FR	15 FR Endo, attached Trocar, J Connector
2015-9619J	19FR	19 FR Endo, attached Trocar, J Connector
2015-9621J	21FR	21 FR Endo, attached Trocar, J Connector




2015-96xxJ

SILICONE ENDO DRAIN (TROCAR/BENDABLE TROCAR + J CONNECTOR)

Endo Drain with Bendable Trocar:


2015-9610B	10FR	10 FR Endo, attached Bendable Trocar
2015-9615B	15FR	15 FR Endo, attached Bendable Trocar
2015-9619B	19FR	19 FR Endo, attached Bendable Trocar
2015-9621B	21FR	21 FR Endo, attached Bendable Trocar



2015-96xxB

Endo Drain Kit (Bendable Trocar + J Connector):

2015-9610JB	10FR	10 FR Endo, attached Bendable Trocar, J Connector
2015-9615JB	15FR	15 FR Endo, attached Bendable Trocar, J Connector
2015-9619JB	19FR	19 FR Endo, attached Bendable Trocar, J Connector
2015-9621JB	21FR	21 FR Endo, attached Bendable Trocar, J Connector




2015-96xxJB

Bendable Trocar:

- Made of stainless steel.
- Bendable, sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

2015-9110	10FR	
2015-9116	16FR	- Bendable
		- Stainless steel
2015-9119	19FR	- Protection cover
2015-9121	21FR	



2015-91xx

SILICONE HUBLESS ENDO DRAIN

Fortune® Silicone Hubless Endo Drain Features:

- Medical grade silicone for superior biocompatibility.
- Hubless design
- X-ray opaque line.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- A great variety of sizes is available.
- Sterilization by E.O. gas.
- Various combination is available.


- A** Hubless Endo Drain
- B** Hubless Endo Drain Kit
- C** Hubless Endo Drain Set*
- D** Hubless Endo Drain System*
- E** Hubless Endo Drain Bag System*



SPECIFICATIONS

Hubless Endo Drain:


REF. No.	Size	Description
2015-0710	10 FR	
2015-0715	15 FR	
2015-0719	19 FR	- Transparent - 1200 mm
2015-0721	21 FR	- X-ray opaque line
2015-0724	24 FR	
2015-0729	29 FR	



2015-07XX

Hubless Endo Drain Kit:

REF. No.	Size	Description
2015-9710	10 FR	10 FR Hubless Endo, attached Trocar
2015-9715	15 FR	15 FR Hubless Endo, attached Trocar
2015-9719	19 FR	19 FR Hubless Endo, attached Trocar
2015-9721	21 FR	21 FR Hubless Endo, attached Trocar




2015-97XX

SILICONE HUBLESS ENDO DRAIN


Hubless Endo Drain Set*:

REF. No.	Size	Description
0610-0710	10 FR	10 FR Hubless Endo, 100 cc Reservoir
0610-0715	15 FR	15 FR Hubless Endo, 100 cc Reservoir
0610-0719	19 FR	19 FR Hubless Endo, 100 cc Reservoir
0610-0721	21 FR	21 FR Hubless Endo, 100 cc Reservoir



Hubless Endo Drain System*:

REF. No.	Size	Description
0610-9710	10 FR	10 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9715	15 FR	15 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9719	19 FR	19 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9721	21 FR	21 FR Hubless Endo, Trocar, 100 cc Reservoir




Hubless Endo Drain Bag System*:

REF. No.	Size	Description
0910-971020	10 FR	10 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-971520	15 FR	15 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-971920	19 FR	19 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-972120	21 FR	21 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag

* Hubless Endo Drain Set, Hubless Endo Drain System and Hubless Endo Drain Bag System are available with 150 cc, 200 cc or 400 cc Reservoir.
For detailed model numbers, please see the attached appendix.

* Hubless Endo Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



SILICONE HUBLESS ENDO DRAIN (TROCAR/BENDABLE TROCAR+ J CONNECTOR)

Trocar Features: Please refer to the introduction on P5-4.

J Connector Features:


J Connector is designed to connect Endo Drain with the size of 10FR/15FR/19FR/21FR and Innerspring Reservoir. It is made of PVC which is a Non-Toxic material, with a little bendable character to reduce the risk of fracture. The side of J Connector designed to connect with Endo Drain is tapering, making assembly and disassembly more easily.



SPECIFICATIONS

Hubless Endo drain + J Connector:


2015-0710J	10FR	10 FR Hubless Endo, attached J Connector
2015-0715J	15FR	15 FR Hubless Endo, attached J Connector
2015-0719J	19FR	19 FR Hubless Endo, attached J Connector
2015-0721J	21FR	21 FR Hubless Endo, attached J Connector



2015-07xxJ

Hubless Endo Drain Kit (Trocar + J Connector):

2015-9710J	10FR	10 FR Hubless Endo, attached Trocar, J Connector
2015-9715J	15FR	15 FR Hubless Endo, attached Trocar, J Connector
2015-9719J	19FR	19 FR Hubless Endo, attached Trocar, J Connector
2015-9721J	21FR	21 FR Hubless Endo, attached Trocar, J Connector




2015-97xxJ

GENERAL SURGICAL SERIES

SILICONE HUBLESS ENDO DRAIN (TROCAR/BENDABLE TROCAR+ J CONNECTOR)

Hubless Endo Drain with Bendable Trocar:


2015-9710B	10FR	10 FR Hubless Endo, attached Bendable Trocar
2015-9715B	15FR	15 FR Hubless Endo, attached Bendable Trocar
2015-9719B	19FR	19 FR Hubless Endo, attached Bendable Trocar
2015-9721B	21FR	21 FR Hubless Endo, attached Bendable Trocar



2015-97xxB

Hubless Endo Drain Kit (Bendable Trocar + J Connector):

2015-9710JB	10FR	10 FR Hubless Endo, attached Bendable Trocar, J Connector
2015-9715JB	15FR	15 FR Hubless Endo, attached Bendable Trocar, J Connector
2015-9719JB	19FR	19 FR Hubless Endo, attached Bendable Trocar, J Connector
2015-9721JB	21FR	21 FR Hubless Endo, attached Bendable Trocar, J Connector



2015-97xxJB

Bendable Trocar:

- Made of stainless steel.
- Bendable, sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

2015-9110	10FR	
2015-9116	16FR	- Bendable
2015-9119	19FR	- Stainless steel
2015-9121	21FR	- Protection cover




2015-91xx

ENDO CONNECTOR


Fortune® Endo Connector Features:

- Endo Connector is intended to connect Endo Drain or Hubless Endo Drain, introducing body fluid drainage from body cavities.




SPECIFICATIONS

REF. No.	Description
2015-8601	1:1 with plugs, connects 1 to 1 Endo-Drain of 19, 21, 24 & 29 FR, source port fits 10-13 mm tubing
2015-8603	3:1 with plugs, connects 1 to 3 Endo-Drain of 19, 21, 24 & 29 FR, source port fits 10-13 mm tubing



2015-8601



2015-8603

**INTENDED USE:**

Fortune wound drainage tubes are intended to drain fluids from body cavities.


Fortune FixPad® Anchoring Device can be used along with different types of “Fortune” wound drainage tubes.

Fortune® Wound Drainage Tubes Features:

- Medical grade silicone for superior biocompatibility.
- Excellent tear and tensile strength, great elongation, flexibility, and durometer range.
- Different drain types are available (Flat Drain/ Vac Drain/ Round Drain/ Endo Drain/ Hubless Endo drain)

FixPad® Anchoring Device Features:


- Applicable for products, i.e. round tube and catheter, the size scope between 3 mm to 12 mm (9 FR to 36 FR).
- The adhesive pad is made of nonwoven fabric and coated with medical pressure sensitive adhesive.
- Easy to apply and remove without sticky residue.
- keep the skin dry and clean before use.
- Less likely to cause skin allergies.




SPECIFICATIONS

Silicone Flat Drain & FixPad® Anchoring Device:

REF. No.	Size	Flat Drain Description
FD14-0051	2 x 4.5 mm	
FD14-0071	3 x 7 mm	3/4 perforation on 200 mm drain X-ray opaque effective
FD14-0101	4 x 10 mm	
FD14-0131	5 x 13 mm	
FD14-0052	2 x 4.5 mm	
FD14-0072	3 x 7 mm	Full perforation on 200 mm drain X-ray opaque effective
FD14-0102	4 x 10 mm	
FD14-0132	5 x 13 mm	



FD14-OXX1 3/4 Perforation





FD14-OXX2 Full Perforation


“FORTUNE” WOUND DRAINAGE TUBE WITH FIXPAD®

**Silicone Vac Drain & FixPad® Anchoring Device:**


REF. No.	Size	Vac Drain Description
FD14-0410	10 FR	
FD14-0415	15 FR	- Full round Vac Drain
FD14-0419	19 FR	- 200 mm round duct drain, 600 mm round tube
FD14-0421	21 FR	- X-ray opaque effective
FD14-0421	21 FR	
FD14-0408	3,5 x 8 mm	- Full flat Vac Drain
FD14-0408	3,5 x 8 mm	- 200 mm flat duct drain, 600 mm round tube
FD14-0411	4 x 11 mm	- X-ray opaque effective

**Silicone Round Drain & FixPad® Anchoring Device:**


REF. No.	Size	Round Drain Description
FD15-0210	10 FR	
FD15-0212	12 FR	
FD15-0214	14 FR	- End perforation
FD15-0216	16 FR	- 900 mm,
FD15-0218	18 FR	- X-ray opaque line
FD15-0220	20 FR	

**Silicone Endo Drain & FixPad® Anchoring Device:**

REF. No.	Size	Endo Drain Description
FD15-0610	10 FR	
FD15-0615	15 FR	
FD15-0619	19 FR	- Transparent
FD15-0621	21 FR	- 1200 mm
FD15-0624	24 FR	- X-ray opaque line
FD15-0629	29 FR	

**Silicone Hubless Endo Drain & FixPad® Anchoring Device:**

REF. No.	Size	Hubless Endo Drain Description
FD15-0710	10 FR	
FD15-0715	15 FR	
FD15-0719	19 FR	- Transparent
FD15-0721	21 FR	- 1200 mm
FD15-0724	24 FR	- X-ray opaque line
FD15-0729	29 FR	



Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Flat Drain Set (Top outlet Reservoir + Flat Drain)	0610-0051	100 cc Reservoir + Flat Drain	10	8	80	3.2	4.5 mm Flat 3/4 + 100 cc Reservoir
	0610-0071						7 mm Flat 3/4 + 100 cc Reservoir
	0610-0101						10 mm Flat 3/4 + 100 cc Reservoir
	0610-0131						13 mm Flat 3/4 + 100 cc Reservoir
	0610-0052						4.5 mm Flat full + 100 cc Reservoir
	0610-0072						7 mm Flat full + 100 cc Reservoir
	0610-0102						10 mm Flat full + 100 cc Reservoir
	0610-0132						13 mm Flat full + 100 cc Reservoir
	0615-0051	150 cc Reservoir + Flat Drain	10	8	80	3.2	4.5 mm Flat 3/4 + 150 cc Reservoir
	0615-0071						7 mm Flat 3/4 + 150 cc Reservoir
	0615-0101						10 mm Flat 3/4 + 150 cc Reservoir
	0615-0131						13 mm Flat 3/4 + 150 cc Reservoir
	0615-0052						4.5 mm Flat full + 150 cc Reservoir
	0615-0072						7 mm Flat full + 150 cc Reservoir
	0615-0102						10 mm Flat full + 150 cc Reservoir
	0615-0132						13 mm Flat full + 150 cc Reservoir
	0620-0051	200 cc Reservoir + Flat Drain	5	8	40	3.2	4.5 mm Flat 3/4 + 200 cc Reservoir
	0620-0071						7 mm Flat 3/4 + 200 cc Reservoir
	0620-0101						10 mm Flat 3/4 + 200 cc Reservoir
	0620-0131						13 mm Flat 3/4 + 200 cc Reservoir
	0620-0052						4.5 mm Flat full + 200 cc Reservoir
	0620-0072						7 mm Flat full + 200 cc Reservoir
	0620-0102						10 mm Flat full + 200 cc Reservoir
	0620-0132						13 mm Flat full + 200 cc Reservoir
	0640-0071	400 cc Reservoir + Flat Drain	5	8	40	4	7 mm Flat 3/4 + 400 cc Reservoir
	0640-0101						10 mm Flat 3/4 + 400 cc Reservoir
	0640-0131						13 mm Flat 3/4 + 400 cc Reservoir
	0640-0072						7 mm Flat full + 400 cc Reservoir
	0640-0102						10 mm Flat full + 400 cc Reservoir
	0640-0132						13 mm Flat full + 400 cc Reservoir
Flat Drain System (Top outlet Reservoir + Flat Drain + Trocar)	0610-9051	100 cc Reservoir + Flat Drain Kit	10	8	80	3.2	4.5 mm Flat 3/4 + Trocar + 100 cc Reservoir
	0610-9071						7 mm Flat 3/4 + Trocar + 100 cc Reservoir
	0610-9101						10 mm Flat 3/4 + Trocar + 100 cc Reservoir
	0610-9131						13 mm Flat 3/4 + Trocar + 100 cc Reservoir
	0610-9052						4.5 mm Flat full + Trocar + 100 cc Reservoir
	0610-9072						7 mm Flat full + Trocar + 100 cc Reservoir
	0610-9102						10 mm Flat full + Trocar + 100 cc Reservoir
	0610-9132						13 mm Flat full + Trocar + 100 cc Reservoir
	0615-9051	150 cc Reservoir + Flat Drain Kit	10	8	80	3.2	4.5 mm Flat 3/4 + Trocar + 150 cc Reservoir
	0615-9071						7 mm Flat 3/4 + Trocar + 150 cc Reservoir
	0615-9101						10 mm Flat 3/4 + Trocar + 150 cc Reservoir
	0615-9131						13 mm Flat 3/4 + Trocar + 150 cc Reservoir
	0615-9052						4.5 mm Flat full + Trocar + 150 cc Reservoir
	0615-9072						7 mm Flat full + Trocar + 150 cc Reservoir
	0615-9102						10 mm Flat full + Trocar + 150 cc Reservoir
	0615-9132						13 mm Flat full + Trocar + 150 cc Reservoir
	0620-9051	200 cc Reservoir + Flat Drain Kit	5	8	40	3.2	4.5 mm Flat 3/4 + Trocar + 200 cc Reservoir
	0620-9071						7 mm Flat 3/4 + Trocar + 200 cc Reservoir
	0620-9101						10 mm Flat 3/4 + Trocar + 200 cc Reservoir
	0620-9131						13 mm Flat 3/4 + Trocar + 200 cc Reservoir
	0620-9052						4.5 mm Flat full + Trocar + 200 cc Reservoir
	0620-9072						7 mm Flat full + Trocar + 200 cc Reservoir
	0620-9102						10 mm Flat full + Trocar + 200 cc Reservoir
	0620-9132						13 mm Flat full + Trocar + 200 cc Reservoir
	0640-9071	400 cc Reservoir + Flat Drain Kit	5	8	40	4	7 mm Flat 3/4 + Trocar + 400 cc Reservoir
	0640-9101						10 mm Flat 3/4 + Trocar + 400 cc Reservoir
	0640-9131						13 mm Flat 3/4 + Trocar + 400 cc Reservoir
	0640-9072						7 mm Flat full + Trocar + 400 cc Reservoir
	0640-9102						10 mm Flat full + Trocar + 400 cc Reservoir
	0640-9132						13 mm Flat full + Trocar + 400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Flat Drain Bag System (Bottom outlet Reservoir + Flat Drain + Trocar + Collection Bag)	0910-905120	100 cc Reservoir + Flat Drain Kit + 200 cc Collection Bag	10	8	80	4	4.5 mm Flat 3/4 + Trocar + 100 Reservoir +200 Bag
	0910-907120						7 mm Flat 3/4 + Trocar + 100 Reservoir + 200 Bag
	0910-910120						10 mm Flat 3/4 + Trocar + 100 Reservoir + 200 Bag
	0910-913120						13 mm Flat 3/4 + Trocar + 100 Reservoir + 200 Bag
	0910-905220						4.5 mm Flat full + Trocar + 100 Reservoir + 200 Bag
	0910-907220						7 mm Flat full + Trocar + 100 Reservoir + 200 Bag
	0910-910220						10 mm Flat full + Trocar + 100 Reservoir + 200 Bag
	0910-913220						13 mm Flat full + Trocar + 100 Reservoir + 200 Bag
	0915-905140	150 cc Reservoir + Flat Drain Kit + 400 cc Collection Bag	10	8	80	4	4.5 mm Flat 3/4 + Trocar + 150 Reservoir + 400 Bag
	0915-907140						7 mm Flat 3/4 + Trocar + 150 Reservoir + 400 Bag
	0915-910140						10 mm Flat 3/4 + Trocar + 150 Reservoir + 400 Bag
	0915-913140						13 mm Flat 3/4 + Trocar + 150 Reservoir + 400 Bag
	0915-905240						4.5 mm Flat full + Trocar + 150 Reservoir + 400 Bag
	0915-907240						7 mm Flat full + Trocar + 150 Reservoir + 400 Bag
	0915-910240						10 mm Flat full + Trocar + 150 Reservoir + 400 Bag
	0915-913240						13 mm Flat full + Trocar + 150 Reservoir + 400 Bag
	0920-905140	200 cc Reservoir + Flat Drain Kit + 400 cc Collection Bag	5	8	40	4	4.5 mm Flat 3/4 + Trocar + 200 Reservoir + 400 Bag
	0920-907140						7 mm Flat 3/4 + Trocar + 200 Reservoir + 400 Bag
	0920-910140						10 mm Flat 3/4 + Trocar + 200 Reservoir + 400 Bag
	0920-913140						13 mm Flat 3/4 + Trocar + 200 Reservoir + 400 Bag
	0920-905240						4.5 mm Flat full + Trocar + 200 Reservoir + 400 Bag
	0920-907240						7 mm Flat full + Trocar + 200 Reservoir + 400 Bag
	0920-910240						10 mm Flat full + Trocar + 200 Reservoir + 400 Bag
	0920-913240						13 mm Flat full + Trocar + 200 Reservoir + 400 Bag
	0940-907180	400 cc Reservoir + Flat Drain Kit + 800 cc Collection Bag	5	8	40	4	7 mm Flat 3/4 + Trocar + 400 Reservoir + 800 Bag
	0940-910180						10 mm Flat 3/4 + Trocar + 400 Reservoir + 800 Bag
	0940-913180						13 mm Flat 3/4 + Trocar + 400 Reservoir + 800 Bag
	0940-907280						7 mm Flat full + Trocar + 400 Reservoir + 800 Bag
	0940-910280						10 mm Flat full + Trocar + 400 Reservoir + 800 Bag
	0940-913280						13 mm Flat full + Trocar + 400 Reservoir + 800 Bag
Vac-Drain Set (Top outlet Reservoir + Vac-Drain)	0610-0410	100 cc Reservoir + full Vac-Drain	10	8	80	3.2	10 FR full Vac + 100 cc Reservoir
	0610-0415						15 FR full Vac + 100 cc Reservoir
	0610-0419						19 FR full Vac + 100 cc Reservoir
	0610-0421						21 FR full Vac + 100 cc Reservoir
	0610-0408						8 mm full Vac + 100 cc Reservoir
	0610-0411	150 cc Reservoir + full Vac-Drain	10	8	80	3.2	11 mm full Vac + 100 cc Reservoir
	0615-0410						10 FR full Vac + 150 cc Reservoir
	0615-0415						15 FR full Vac + 150 cc Reservoir
	0615-0419						19 FR full Vac + 150 cc Reservoir
	0615-0421						21 FR full Vac + 150 cc Reservoir
	0615-0408						8 mm full Vac + 150 cc Reservoir
	0615-0411	200 cc Reservoir + full Vac-Drain	5	8	40	3.2	11 mm full Vac + 150 cc Reservoir
	0620-0410						10 FR full Vac + 200 cc Reservoir
	0620-0415						15 FR full Vac + 200 cc Reservoir
	0620-0419						19 FR full Vac + 200 cc Reservoir
	0620-0421						21 FR full Vac + 200 cc Reservoir
	0620-0408						8 mm full Vac + 200 cc Reservoir
	0620-0411	400 cc Reservoir + full Vac-Drain	5	8	40	4	11 mm full Vac + 200 cc Reservoir
	0640-0415						15 FR full Vac + 400 cc Reservoir
	0640-0419						19 FR full Vac + 400 cc Reservoir
	0640-0421						21 FR full Vac + 400 cc Reservoir
	0640-0408						8 mm full Vac + 400 cc Reservoir
	0640-0411						11 mm full Vac + 400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Vac-Drain Sytem (Top outlet Reservoir + Vac-Drain + Trocar)	0610-9410 0610-9415 0610-9419 0610-9421	100 cc Reservoir + full Vac-Drain Kit	10	8	80	3.2	10 FR full Vac + Trocar + 100 cc Reservoir 15 FR full Vac + Trocar + 100 cc Reservoir 19 FR full Vac + Trocar + 100 cc Reservoir 21 FR full Vac + Trocar + 100 cc Reservoir
	0610-9408 0610-9411						8 mm full Vac + Trocar + 100 cc Reservoir 11 mm full Vac + Trocar + 100 cc Reservoir
	0615-9410 0615-9415 0615-9419 0615-9421	150 cc Reservoir + full Vac-Drain Kit	10	8	80	3.2	10 FR full Vac + Trocar + 150 cc Reservoir 15 FR full Vac + Trocar + 150 cc Reservoir 19 FR full Vac + Trocar + 150 cc Reservoir 21 FR full Vac + Trocar + 150 cc Reservoir
	0615-9408 0615-9411						8 mm full Vac + Trocar + 150 cc Reservoir 11 mm full Vac + Trocar + 150 cc Reservoir
	0620-9410 0620-9415 0620-9419 0620-9421	200 cc Reservoir + full Vac-Drain Kit	5	8	40	4	10 FR full Vac + Trocar + 200 cc Reservoir 15 FR full Vac + Trocar + 200 cc Reservoir 19 FR full Vac + Trocar + 200 cc Reservoir 21 FR full Vac + Trocar + 200 cc Reservoir
	0620-9408 0620-9411						8 mm full Vac + Trocar + 200 cc Reservoir 11 mm full Vac + Trocar + 200 cc Reservoir
	0640-9415 0640-9419 0640-9421	400 cc Reservoir + full Vac-Drain Kit	5	8	40	4	15 FR full Vac + Trocar + 400 cc Reservoir 19 FR full Vac + Trocar + 400 cc Reservoir 21 FR full Vac + Trocar + 400 cc Reservoir
	0640-9408 0640-9411						8 mm full Vac + Trocar + 400 Reservoir 11 mm full Vac + Trocar+ 400 cc Reservoir
Vac-Drain Bag Sytem (Bottom outlet Reservoir + Vac-Drain + Trocar + Collection Bag)	0910-941020 0910-941520 0910-941920 0910-942120	100 cc Reservoir + full Vac-Drain Kit + 200 cc Collection Bag	10	8	80	4	10 FR full Vac + Trocar + 100 Reservoir + 200 Bag 15 FR full Vac + Trocar + 100 Reservoir + 200 Bag 19 FR full Vac + Trocar + 100 Reservoir + 200 Bag 21 FR full Vac + Trocar + 100 Reservoir + 200 Bag
	0910-940820 0910-941120						8 mm full Vac + Trocar + 100 Reservoir + 200 Bag 11 mm full Vac + Trocar + 100 Reservoir + 200 Bag
	0915-941040 0915-941540 0915-941940 0915-942140	150 cc Reservoir + full Vac-Drain Kit + 400 cc Collection Bag	10	8	80	4	10 FR full Vac + Trocar + 150 Reservoir + 400 Bag 15 FR full Vac + Trocar + 150 Reservoir + 400 Bag 19 FR full Vac + Trocar + 150 Reservoir + 400 Bag 21 FR full Vac + Trocar + 150 Reservoir + 400 Bag
	0915-940840 0915-941140						8 mm full Vac + Trocar + 150 Reservoir + 400 Bag 11 mm full Vac + Trocar + 150 Reservoir + 400 Bag
	0920-941040 0920-941540 0920-941940 0920-942140	200 cc Reservoir + full Vac-Drain Kit + 400 cc Collection Bag	5	8	40	4	10 FR full Vac + Trocar + 200 Reservoir + 400 Bag 15 FR full Vac + Trocar + 200 Reservoir + 400 Bag 19 FR full Vac + Trocar + 200 Reservoir + 400 Bag 21 FR full Vac + Trocar + 200 Reservoir + 400 Bag
	0920-940840 0920-941140						8 mm full Vac + Trocar + 200 Reservoir + 400 Bag 11 mm full Vac + Trocar + 200 Reservoir + 400 Bag
	0940-941580 0940-941980 0940-942180	400 cc Reservoir + full Vac-Drain Kit + 800 cc Collection Bag	5	8	40	4	15 FR full Vac + Trocar + 400 Reservoir + 800 Bag 19 FR full Vac + Trocar + 400 Reservoir + 800 Bag 21 FR full Vac + Trocar + 400 Reservoir + 800 Bag
	0940-940880 0940-941180						8 mm full Vac + Trocar + 400 Reservoir + 800 Bag 11 mm full Vac + Trocar + 400 Reservoir + 800 Bag
Silicone Endo-Drain Set (Top outlet Reservoir + Endo-Drain)	0610-0610 0610-0615	100 cc Reservoir + Endo-Drain	10	8	80	3.2	10 FR Endo + 100 cc Reservoir 15 FR Endo + 100 cc Reservoir
	0610-0619 0610-0621						19 FR Endo + 100 cc Reservoir 21 FR Endo + 100 cc Reservoir
	0615-0610 0615-0615 0615-0619 0615-0621	150 cc Reservoir + Endo-Drain	10	8	80	3.2	10 FR Endo + 150 cc Reservoir 15 FR Endo + 150 cc Reservoir 19 FR Endo + 150 cc Reservoir 21 FR Endo + 150 cc Reservoir
	0620-0610 0620-0615 0620-0619 0620-0621						10 FR Endo + 200 cc Reservoir 15 FR Endo + 200 cc Reservoir 19 FR Endo + 200 cc Reservoir 21 FR Endo + 200 cc Reservoir
	0640-0615 0640-0619 0640-0621	400 cc Reservoir + Endo-Drain	5	8	40	4	15 FR Endo + 400 cc Reservoir 19 FR Endo + 400 cc Reservoir 21 FR Endo + 400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Silicone Endo-Drain System (Top outlet Reservoir + Endo-Drain + Trocar)	0610-9610	100 cc Reservoir + Endo-Drain Kit	10	8	80	3.2	10 FR Endo + Trocar + 100 cc Reservoir
	0610-9615						15 FR Endo + Trocar + 100 cc Reservoir
	0610-9619						19 FR Endo + Trocar + 100 cc Reservoir
	0610-9621						21 FR Endo + Trocar + 100 cc Reservoir
	0615-9610	150 cc Reservoir + Endo-Drain Kit	10	8	80	3.2	10 FR Endo + Trocar + 150 cc Reservoir
	0615-9615						15 FR Endo + Trocar + 150 cc Reservoir
	0615-9619						19 FR Endo + Trocar + 150 cc Reservoir
	0615-9621						21 FR Endo + Trocar + 150 cc Reservoir
	0620-9610	200 cc Reservoir + Endo-Drain Kit	5	8	40	4	10 FR Endo + Trocar + 200 cc Reservoir
	0620-9615						15 FR Endo + Trocar + 200 cc Reservoir
	0620-9619						19 FR Endo + Trocar + 200 cc Reservoir
	0620-9621						21 FR Endo + Trocar + 200 cc Reservoir
	0640-9615	400 cc Reservoir + Endo-Drain Kit	5	8	40	4	15 FR Endo + Trocar + 400 cc Reservoir
	0640-9619						19 FR Endo + Trocar + 400 cc Reservoir
	0640-9621						21 FR Endo + Trocar + 400 cc Reservoir
Silicone Endo-Drain Bag System (Bottom outlet Reservoir + Endo-Drain + Trocar + Collection Bag)	0910-961020	100 cc Reservoir + Endo-Drain Kit + 200 cc Bag	10	8	80	4	10 FR Endo + Trocar + 100 Reservoir + 200 Bag
	0910-961520						15 FR Endo + Trocar + 100 Reservoir + 200 Bag
	0910-961920						19 FR Endo + Trocar + 100 Reservoir + 200 Bag
	0910-962120						21 FR Endo + Trocar + 100 Reservoir + 200 Bag
	0915-961040	150 cc Reservoir + Endo-Drain Kit+ 400 cc Bag	10	8	80	4	10 FR Endo + Trocar + 150 Reservoir + 400 Bag
	0915-961540						15 FR Endo + Trocar + 150 Reservoir + 400 Bag
	0915-961940						19 FR Endo + Trocar + 150 Reservoir + 400 Bag
	0915-962140						21 FR Endo + Trocar + 150 Reservoir + 400 Bag
	0920-961040	200 cc Reservoir + Endo-Drain Kit + 400 cc Bag	5	8	40	4	10 FR Endo + Trocar + 200 Reservoir + 400 Bag
	0920-961540						15 FR Endo + Trocar + 200 Reservoir + 400 Bag
	0920-961940						19 FR Endo + Trocar + 200 Reservoir + 400 Bag
	0920-962140						21 FR Endo + Trocar + 200 Reservoir + 400 Bag
	0940-961580	400 cc Reservoir + Endo-Drain Kit + 800 cc Bag	5	8	40	4	15 FR Endo + Trocar + 400 Reservoir + 800 Bag
	0940-961980						19 FR Endo + Trocar + 400 Reservoir + 800 Bag
	0940-962180						21 FR Endo + Trocar + 400 Reservoir + 800 Bag
Silicone Hubless Endo-Drain Set (Top outlet Reservoir + Endo-Drain)	0610-0710	100 cc Reservoir + Endo-Drain	10	8	80	3.2	10 FR Hubless Endo + 100 cc Reservoir
	0610-0715						15 FR Hubless Endo + 100 cc Reservoir
	0610-0719						19 FR Hubless Endo + 100 cc Reservoir
	0610-0721						21 FR Hubless Endo + 100 cc Reservoir
	0615-0710	150 cc Reservoir + Endo-Drain	10	8	80	3.2	10 FR Hubless Endo + 150 cc Reservoir
	0615-0715						15 FR Hubless Endo + 150 cc Reservoir
	0615-0719						19 FR Hubless Endo + 150 cc Reservoir
	0615-0721						21 FR Hubless Endo + 150 cc Reservoir
	0620-0710	200 cc Reservoir + Endo-Drain	5	8	40	3.2	10 FR Hubless Endo + 200 cc Reservoir
	0620-0715						15 FR Hubless Endo + 200 cc Reservoir
	0620-0719						19 FR Hubless Endo + 200 cc Reservoir
	0620-0721						21 FR Hubless Endo + 200 cc Reservoir
	0640-0715	400 cc Reservoir + Endo-Drain	5	8	40	4	15 FR Hubless Endo + 400 cc Reservoir
	0640-0719						19 FR Hubless Endo + 400 cc Reservoir
	0640-0721						21 FR Hubless Endo + 400 cc Reservoir
Silicone Hubless Endo-Drain System (Top outlet Reservoir + Endo-Drain + Trocar)	0610-9710	100 cc Reservoir + Endo-Drain Kit	10	8	80	3.2	10 FR Hubless Endo+Trocar+100 cc Reservoir
	0610-9715						15 FR Hubless Endo+Trocar+100 cc Reservoir
	0610-9719						19 FR Hubless Endo+Trocar+100 cc Reservoir
	0610-9721						21 FR Hubless Endo+Trocar+100 cc Reservoir
	0615-9710	150 cc Reservoir + Endo-Drain Kit	10	8	80	3.2	10 FR Hubless Endo+Trocar+150 cc Reservoir
	0615-9715						15 FR Hubless Endo+Trocar+150 cc Reservoir
	0615-9719						19 FR Hubless Endo+Trocar+150 cc Reservoir
	0615-9721						21 FR Hubless Endo+Trocar+150 cc Reservoir
	0620-9710	200 cc Reservoir + Endo-Drain Kit	5	8	40	4	10 FR Hubless Endo+Trocar+200 cc Reservoir
	0620-9715						15 FR Hubless Endo+Trocar+200 cc Reservoir
	0620-9719						19 FR Hubless Endo+Trocar+200 cc Reservoir
	0620-9721						21 FR Hubless Endo+Trocar+200 cc Reservoir
	0640-9715	400 cc Reservoir + Endo-Drain Kit	5	8	40	4	15 FR Hubless Endo+Trocar+400 cc Reservoir
	0640-9719						19 FR Hubless Endo+Trocar+400 cc Reservoir
	0640-9721						21 FR Hubless Endo+Trocar+400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Silicone Hubless Endo- Drain Bag System (Bottom outlet Reservoir + Endo-Drain + Trocar + Collection Bag)	0910-971020	100 cc Reservoir + Endo- Drain Kit + 200 cc Bag	10	8	80	4	10 FR Hubless Endo+Trocar+100 Reservoir+200 Bag
	0910-971520						15 FR Hubless Endo+Trocar+100 Reservoir+200 Bag
	0910-971920						19 FR Hubless Endo+Trocar+100 Reservoir+200 Bag
	0910-972120						21 FR Hubless Endo+Trocar+100 Reservoir+200 Bag
	0915-971040	150 cc Reservoir + Endo- Drain Kit+ 400 cc Bag	10	8	80	4	10 FR Hubless Endo+Trocar+150 Reservoir+400 Bag
	0915-971540						15 FR Hubless Endo+Trocar+150 Reservoir+400 Bag
	0915-971940						19 FR Hubless Endo+Trocar+150 Reservoir+400 Bag
	0915-972140						21 FR Hubless Endo+Trocar+150 Reservoir+400 Bag
	0920-971040	200 cc Reservoir + Endo- Drain Kit + 400 cc Bag	5	8	40	4	10 FR Hubless Endo+Trocar+200 Reservoir+400 Bag
	0920-971540						15 FR Hubless Endo+Trocar+200 Reservoir+400 Bag
	0920-971940						19 FR Hubless Endo+Trocar+200 Reservoir+400 Bag
	0920-972140						21 FR Hubless Endo+Trocar+200 Reservoir+400 Bag
	0940-971580	400 cc Reservoir + Endo- Drain Kit + 800 cc Bag	5	8	40	4	15 FR Hubless Endo+Trocar+400 Reservoir+800 Bag
	0940-971980						19 FR Hubless Endo+Trocar+400 Reservoir+800 Bag
	0940-972180						21 FR Hubless Endo+Trocar+400 Reservoir+800 Bag

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Round Drain Set (Top outlet Reservoir + Round Drain)	0610-0208	100 cc Reservoir + Round Drain	10	8	80	3.2	8 FR Round end+ 100 cc Reservoir
	0610-0210						10 FR Round end + 100 cc Reservoir
	0610-0212						12 FR Round end + 100 cc Reservoir
	0610-0214						14 FR Round end + 100 cc Reservoir
	0610-0216						16 FR Round end + 100 cc Reservoir
	0610-0218						18 FR Round end + 100 cc Reservoir
	0610-0220						20 FR Round end + 100 cc Reservoir
	0610-0308						8 FR Round mid + 100 cc Reservoir
	0610-0310						10 FR Round mid + 100 cc Reservoir
	0610-0312						12 FR Round mid + 100 cc Reservoir
	0610-0314						14 FR Round mid + 100 cc Reservoir
	0610-0316						16 FR Round mid + 100 cc Reservoir
	0610-0318						18 FR Round mid + 100 cc Reservoir
	0610-0320						20 FR Round mid + 100 cc Reservoir
	0615-0208	150 cc Reservoir + Round Drain	10	8	80	3.2	8 FR Round end + 150 cc Reservoir
	0615-0210						10 FR Round end + 150 cc Reservoir
	0615-0212						12 FR Round end + 150 cc Reservoir
	0615-0214						14 FR Round end + 150 cc Reservoir
	0615-0216						16 FR Round end + 150 cc Reservoir
	0615-0218						18 FR Round end + 150 cc Reservoir
	0615-0220						20 FR Round end + 150 cc Reservoir
	0615-0308						8 FR Round mid + 150 cc Reservoir
	0615-0310						10 FR Round mid + 150 cc Reservoir
	0615-0312						12 FR Round mid + 150 cc Reservoir
	0615-0314						14FR Round mid + 150cc Reservoir
	0615-0316						16 FR Round mid + 150 cc Reservoir
	0615-0318						18 FR Round mid + 150 cc Reservoir
	0615-0320						20 FR Round mid + 150 cc Reservoir
	0620-0208	200 cc Reservoir + Round Drain	5	8	40	3.2	8 FR Round end + 200 cc Reservoir
	0620-0210						10 FR Round end + 200 cc Reservoir
	0620-0212						12 FR Round end + 200 cc Reservoir
	0620-0214						14FR Round end + 200 cc Reservoir
	0620-0216						16 FR Round end + 200 cc Reservoir
	0620-0218						18 FR Round end + 200 cc Reservoir
	0620-0220						20 FR Round end + 200 cc Reservoir
	0620-0308						8 FR Round mid + 200 cc Reservoir
	0620-0310						10 FR Round mid + 200 cc Reservoir
	0620-0312						12 FR Round mid + 200 cc Reservoir
	0620-0314						14 FR Round mid + 200 cc Reservoir
	0620-0316						16 FR Round mid + 200 cc Reservoir
	0620-0318						18 FR Round mid + 200 cc Reservoir
	0620-0320						20 FR Round mid + 200 cc Reservoir
	0640-0216	400 cc Reservoir + Round Drain	5	8	40	4	16 FR Round end + 400 cc Reservoir
	0640-0218						18 FR Round end + 400 cc Reservoir
	0640-0220						20 FR Round end + 400 cc Reservoir
	0640-0316						16 FR Round mid + 400 cc Reservoir
	0640-0318						18 FR Round mid + 400 cc Reservoir
	0640-0320						20 FR Round mid + 400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Round Drain System (Top outlet Reservoir + Round Drain + Trocar)	0610-9208	100 cc Reservoir + Round Drain Kit	10	8	80	3.2	8 FR Round end + Trocar + 100 cc Reservoir
	0610-9210						10 FR Round end + Trocar + 100 cc Reservoir
	0610-9214						14 FR Round end + Trocar + 100 cc Reservoir
	0610-9216						16 FR Round end + Trocar + 100 cc Reservoir
	0610-9218						18 FR Round end + Trocar + 100 cc Reservoir
	0610-9220						20 FR Round end + Trocar + 100 cc Reservoir
	0610-9308						8 FR Round mid + Trocar + 100 cc Reservoir
	0610-9310						10 FR Round mid + Trocar + 100 cc Reservoir
	0610-9314						14 FR Round mid + Trocar + 100 cc Reservoir
	0610-9316						16 FR Round mid + Trocar + 100 cc Reservoir
	0610-9318						18 FR Round mid + Trocar + 100 cc Reservoir
	0610-9320						20 FR Round mid + Trocar + 100 cc Reservoir
	0615-9208	150 cc Reservoir + Round Drain Kit	10	8	80	3.2	8 FR Round end + Trocar + 150 cc Reservoir
	0615-9210						10 FR Round end + Trocar + 150 cc Reservoir
	0615-9214						14 FR Round end + Trocar + 150 cc Reservoir
	0615-9216						16 FR Round end + Trocar + 150 cc Reservoir
	0615-9218						18 FR Round end + Trocar + 150 cc Reservoir
	0615-9220						20 FR Round end + Trocar + 150 cc Reservoir
	0615-9308						8 FR Round mid + Trocar + 150 cc Reservoir
	0615-9310						10 FR Round mid + Trocar + 150 cc Reservoir
	0615-9314						14 FR Round mid + Trocar + 150 cc Reservoir
	0615-9316						16 FR Round mid + Trocar + 150 cc Reservoir
	0615-9318						18 FR Round mid + Trocar + 150 cc Reservoir
	0615-9320						20 FR Round mid + Trocar + 150 cc Reservoir
	0620-9208	200 cc Reservoir + Round Drain Kit	5	8	40	3.2	8 FR Round end + Trocar + 200 cc Reservoir
	0620-9210						10 FR Round end + Trocar + 200 cc Reservoir
	0620-9214						14 FR Round end + Trocar + 200 cc Reservoir
	0620-9216						16 FR Round end + Trocar + 200 cc Reservoir
	0620-9218						18 FR Round end + Trocar + 200 cc Reservoir
	0620-9220						20 FR Round end + Trocar + 200 cc Reservoir
	0620-9308						8 FR Round mid + Trocar + 200 cc Reservoir
	0620-9310						10 FR Round mid + Trocar + 200 cc Reservoir
	0620-9314						14 FR Round mid + Trocar + 200 cc Reservoir
	0620-9316						16 FR Round mid + Trocar + 200 cc Reservoir
	0620-9318						18 FR Round mid + Trocar + 200 cc Reservoir
	0620-9320						20 FR Round mid + Trocar + 200 cc Reservoir
	0640-9216	400 cc Reservoir + Round Drain Kit	5	8	40	4	16 FR Round end + Trocar + 400 cc Reservoir
	0640-9218						18 FR Round end + Trocar + 400 cc Reservoir
	0640-9220						20 FR Round end + Trocar + 400 cc Reservoir
	0640-9316						16 FR Round mid + Trocar + 400 cc Reservoir
	0640-9318						18 FR Round mid + Trocar + 400 cc Reservoir
	0640-9320						20 FR Round mid + Trocar + 400 cc Reservoir

Product	REF. No.	Size	pc /box	box /ctn	pc /ctn	cuft /ctn	Description
Round Drain Bag System (Bottom outlet Reservoir + Round Drain + Trocar+ Collection Bag)	0910-920820	100 cc Reservoir + Round Drain Kit + 200 cc Collection Bag	10	8	80	4	8 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-921020						10 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-921420						14 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-921620						16 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-921820						18 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-922020						20 FR Round end + Trocar + 100 Reservoir + 200 Bag
	0910-930820						8 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0910-931020						10 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0910-931420						14 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0910-931620						16 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0910-931820						18 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0910-932020						20 FR Round mid + Trocar + 100 Reservoir + 200 Bag
	0915-920840	150 cc Reservoir + Round Drain Kit + 400 cc Collection Bag	10	8	80	4	8 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-921040						10 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-921440						14 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-921640						16 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-921840						18 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-922040						20 FR Round end + Trocar + 150 Reservoir + 400 Bag
	0915-930840						8 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0915-931040						10 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0915-931440						14 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0915-931640						16 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0915-931840						18 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0915-932040						20 FR Round mid + Trocar + 150 Reservoir + 400 Bag
	0920-920840		5	8	40	4	8 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-921040						10 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-921440						14 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-921640						16 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-921840						18 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-922040						20 FR Round end + Trocar + 200 Reservoir + 400 Bag
	0920-930840						8 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0920-931040						10 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0920-931440						14 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0920-931640						16 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0920-931840						18 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0920-932040						20 FR Round mid + Trocar + 200 Reservoir + 400 Bag
	0940-921680	400 cc Reservoir + Round Drain Kit + 800 cc Collection Bag	5	8	40	4	16 FR Round end + Trocar + 400 Reservoir + 800 Bag
	0940-921880						18 FR Round end + Trocar + 400 Reservoir + 800 Bag
	0940-922080						20 FR Round end + Trocar + 400 Reservoir + 800 Bag
	0940-931680						16 FR Round mid + Trocar + 400 Reservoir + 800 Bag
	0940-931880						18 FR Round mid + Trocar + 400 Reservoir + 800 Bag
	0940-932080						20 FR Round mid + Trocar + 400 Reservoir + 800 Bag