

SILICONE **VESSEL ID** LOOPS






SILICONE VESSEL ID LOOPS

INTENDED USE:

The "Fortune®" Silicone Vessel ID Loops are used in complicated surgical procedures by offering retraction, occlusion, and identification of arteries, nerves, veins, tendons or ureters.

Fortune® Silicone Vessel ID Loops features:

- Medical grade silicone for superior biocompatibility.
- Soft and smooth material for better comfort
- Silicone radiopaque
- Non-adherent and non-absorbent to tissue
- Non-toxic and latex-free
- Sterile package
- Available in four sizes and four colors





SPECIFICATIONS

REF. No.	Size.	Application	REF. No.	Size.	Application
9005-0503R		Red Color	9005-1309R		Red Color
9005-0503B	Width: 0.5 mm Thickness: 0.3 mm	Blue Color	9005-1309B	Width: 1.3 mm —— Thickness: 0.9 mm	Blue Color
9005-0503Y	Length: 450 mm	Yellow Color	9005-1309Y	Length: 450 mm	Yellow Color
9005-0503W		White Color	9005-1309W		White Color
					<u></u>

9005-0503X 9005-1309X


SILICONE VESSEL ID LOOPS

9005-2510X

REF. No.	Size.	Application	REF. No.	Size.	Application
9005-2510R		Red Color	9005-5013R		Red Color
9005-2510B		Blue Color	9005-5013B	Width: 5.0 mm Thickness: 1.3 mm	Blue Color
9005-2510Y	Thickness: 1.0 mm —— Length: 450 mm	Yellow Color	9005-5013Y	Length: 450 mm	Yellow Color
9005-2510W		White Color	9005-5013W		White Color
9005-2510A		Red, Blue, Yellow, White Color	-		
		} &			



6 F., No. 29, Sec. 2, Jhongjheng E. Rd., Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com



9005-5013X