

SILICONE THORACIC DRAIN

Superior Compatibility for Chest Drainage



THORACIC
SURGICAL
SERIES

SILICONE THORACIC DRAIN


Fortune® Silicone Thoracic Drain Features:

- 100% medical grade silicone for superior biocompatibility.
- Superior strength reduces risk of tearing or breaking.
- Clear, quality silicone provides inspection and drainage visibility.
- Soft, kink-resistant silicone conforms easily to anatomy for improved patient comfort.
- X-ray opaque line.
- Beveled tip and smoothed shaft can minimize trauma during insertion and withdrawal.
- Enlarged drainage holes maximize drainage and prevent clog.
- Expanded/expandable funnel connects easily with larger O.D. connector for further fluid collection.
- Three types are available:
 - S type: Straight type
 - SF type: Straight type with funnel
 - RF type: Right angle type with funnel
- Sterilized double packaging.

Silicone Thoracic Drain, S Type

Silicone Thoracic Drain, SF Type

Silicone Thoracic Drain, RF Type



SPECIFICATIONS

Thoracic Drain, SF Type:

REF. No.	Size	O. D.	I. D.	Description
2018-0116	16 FR	5.3 mm	3.5 mm	
2018-0120	20 FR	6.7 mm	4.5 mm	
2018-0124	24 FR	8.0 mm	5.5 mm	- Straight type with funnel
2018-0128	28 FR	9.3 mm	6.5 mm	- 500 mm
2018-0132	32 FR	10.7 mm	7.5 mm	- X-ray opaque line
2018-0136	36 FR	12.0 mm	8.5 mm	
2018-0140	40 FR	13.3 mm	9.5 mm	




2018-01XX

SILICONE THORACIC DRAIN

Thoracic Drain, S Type:


REF. No.	Size	O. D.	I. D.	Description
2018-0009	9 FR	3.0 mm	1.8 mm	
2018-0012	12 FR	4.0 mm	2.5 mm	
2018-0015	15 FR	5.0 mm	3.3 mm	
2018-0016	16 FR	5.3 mm	3.5 mm	
2018-0018	18 FR	6.0 mm	4.0 mm	
2018-0020	20 FR	6.7 mm	4.5 mm	
2018-0021	21 FR	7.0 mm	4.8 mm	
2018-0022	22 FR	7.3 mm	5.0 mm	- Straight type
2018-0024	24 FR	8.0 mm	5.5 mm	- 500 mm
				- X-ray opaque line
2018-0027	27 FR	9.0 mm	6.3 mm	
2018-0028	28 FR	9.3 mm	6.5 mm	
2018-0030	30 FR	10.0 mm	7.0 mm	
2018-0032	32 FR	10.7 mm	7.5 mm	
2018-0033	33 FR	11.0 mm	7.8 mm	
2018-0034	34 FR	11.3 mm	8.0 mm	
2018-0036	36 FR	12.0 mm	8.5 mm	
2018-0040	40 FR	13.3 mm	9.5 mm	



2018-00XX

Thoracic Drain, RF Type:

REF. No.	Size	O. D.	I. D.	Description
2018-0516	16 FR	5.3 mm	3.5 mm	
2018-0520	20 FR	6.7 mm	4.5 mm	
2018-0524	24 FR	8.0 mm	5.5 mm	- Right angle type with funnel
2018-0528	28 FR	9.3 mm	6.5 mm	- 500 mm
				- X-ray opaque line
2018-0532	32 FR	10.7 mm	7.5 mm	
2018-0536	36 FR	12.0 mm	8.5 mm	
2018-0540	40 FR	13.3 mm	9.5 mm	



2018-05XX



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
 Danshuei Dist, New Taipei City. 251 Taiwan
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
 www.fortunemed.com Mail:info@fortunemed.com